**SANTA  
CASA**  
de Curitiba

# ***Monitoramento da Limpeza***

***Enf<sup>a</sup>. Cristiane Uyeno***

# Responsabilidade

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# *Aumento de 3x*

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Table D.1—In-use tests available to assess efficacy of cleaning of medical devices

Test method	Soil component tested	Limit of detection	Limitations	Length of test (after sample collection)
Ortho-phthaldialdehyde (OPA) method (Fengler, et al., 2001; Verjat, et al., 1999). Swab device or elute device with liquid, then test sample using OPA method.	Protein	0.01 µg/mL	Sensitivity unrealistic (i.e., routine handling with hands could trigger positive reaction).	~ 1 to 5 min
Biuret reaction (Kruger, 1997). Swab device, immerse in reagent, and assess for color development.	Protein	5.5 µg/cm <sup>2</sup>	Not applicable to lumens. Author suggests that > 20 µg/cm <sup>2</sup> is unacceptably high for protein, but no rationale is given for this benchmark. Rust causes color interference.	10 min
Protein method. Swab device, immerse in reagent, and assess for color development.	Protein	Not indicated	Not applicable to lumens. No indication is given of what level of soiling is present for a positive test result.	Stated as "minutes"
ATP method. Swab device, extract ATP from swab, and determine ATP. Or use fluid rinse as sample.	ATP (present in eukaryotic cells and live bacteria)	Not indicated	Needs instrumentation to read test. Requires cells (eukaryotic or prokaryotic) to be present. No ATP is detected if only protein or carbohydrate is present.	30 seconds
ATP bioluminescence (Davidson, et al., 1999).	Bacteria ( <i>S. aureus</i> and <i>E. coli</i> )	<10 <sup>4</sup> cfu/100 cm <sup>2</sup>	Not indicated.	Stated as "minutes"
Ninhydrin test (deBruin, 2002). Swab device, immerse swab in test reagent, and assess for color development.	Protein	2.5 µg/swab	Not applicable to lumens; interference in color detection by rust, etc., from cleaned devices that mask swab color.	20 min
UV-VIS spectroscopy (Kneiler, 2001).	Residual blood	Not indicated	Not indicated.	Not indicated
Limulus amoebocyte lysate assay (LAL). Elute device with liquid, then test sample using LAL method.	Endotoxin	0.0032 EU/mL	Sensitivity unrealistic (i.e., routine handling could trigger positive reaction); does not detect proteins, organic matter, or viable microorganisms.	10 to 30 min
Expose cleaned medical device to a 2 % hydrogen peroxide solution. If the solution bubbles, then there is residual material on the device that contains catalase. Because blood or other cellular components are the most likely source, the presence of bubbles indicates that cleaning was inadequate. (AAMI TIR12)	Catalase containing material (e.g., eukaryotic cells and some bacterial cells)	Not indicated	Instruments must be further cleaned to remove hydrogen peroxide residue.	~ 1 to 5 min

NOTE—More comprehensive lists of test soils and test methods are provided in Tables 5 and 6 of AAMI TIR30.

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




# Testes

If the liquid or swab is green, the surface is "clean", i.e. the swab contains less than 3µg of protein residues.

If the liquid or swab is any shade of grey or purple, protein residue is present. The surface should be re-cleaned and re-tested.

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